

# BIODEX ADVANTAGE SOFTWARE FOR MULTI-JOINT SYSTEM 4

---

VERSION 4.56

830-000  
835-000  
840-000  
850-000



## BIODEX


Biomedex Medical Systems, Inc.

20 Ramsey Road, Shirley, New York, 11967-4704, Tel: 800-224-6339 (toll free 631-924-9000), Fax: 631-924-9241, Email: info@biomedex.com, www.biomedex.com

---

The forthcoming System 4 Advantage software release 4.56 will include the following:

1. The new enhanced data export utility provides a simple means to export patient data files of multiple subjects' calculated data and raw data from tests in the Advantage database all at once. Calculated data can be exported as .csv or .txt file formats. "Raw" data files are exported as .txt files which are compatible with popular third party research programs, such as MatLab, LabView or BioPac.



2. This utility is available in the Biodex Medical program folder and can be launched by clicking on the Multi-Joint Export icon.
3. The Rehab Exported Data Parser, in the form of an Excel macro, facilitates the process of working with exported data. The Data parser has the header files preformatted. The data can also be presented in a similar format as a report.
4. Expanded and updated Isokinetic normative data is available for knee, ankle, hip, shoulder, elbow, and wrist. Peak Torque to BW% for each direction (away and towards) as well as a Goal for agonist/antagonist ratio (which is the normally weaker muscle group/normally stronger group) from a population of 178 (93 male and 85 female) healthy non-athletic volunteers aged 15–83 years\*. The Biodex Multi-joint System now includes normative data from pediatric age (6 -13) years through adults to age 83.

## ISOKINETIC NORMATIVE DATA

Gravity Correction: Yes    Windowing: Yes

### KNEE EXTENSION/FLEXION – 90 degrees per second

Agonist/Antagonist Ratio = Normally Weaker/Normally Stronger

*Peak Torque to Body Weight %*

**AWAY**

**Extension**

Gender : Age	<30	30-39	40-49	50-59	60-69	>70
male high	116%	96%	91%	90%	82%	71%
male low	79%	78%	67%	57%	56%	55%
female high	88%	88%	79%	71%	66%	56%
female low	69%	64%	61%	53%	42%	43%

**TOWARD**

**Flexion**

Gender : Age	<30	30-39	40-49	50-59	60-69	>70
male high	62%	47%	46%	59%	47%	38%
male low	35%	32%	32%	28%	29%	24%
female high	48%	42%	39%	38%	33%	31%
female low	38%	29%	25%	25%	20%	21%

**Agonist/Antagonist Ratio: Flexion/Extension Ratio**

Gender : Age	<30	30-39	40-49	50-59	60-69	>70
male high	53%	49%	50%	65%	57%	54%
male low	44%	41%	47%	50%	52%	44%
average	48%	45%	49%	57%	54%	49%
female high	55%	47%	49%	54%	50%	55%
female low	55%	45%	42%	47%	47%	49%
average	55%	46%	45%	50%	48%	52%

### ANKLE PLANTAR/DORSIFLEXION – 60 degrees per second

*Peak Torque to Body Weight %*

**AWAY**

**Plantarflexion**

Gender : Age	<30	30-39	40-49	50-59	60-69	>70
male high	69%	57%	55%	57%	53%	46%
male low	47%	39%	37%	41%	34%	25%
female high	90%	61%	48%	68%	61%	45%
female low	47%	48%	45%	41%	36%	31%

**TOWARD**

**Dorsiflexion**

Gender : Age	<30	30-39	40-49	50-59	60-69	>70
male high	19%	16%	18%	17%	15%	13%
male low	13%	15%	14%	15%	14%	12%
female high	14%	15%	13%	14%	12%	11%
female low	11%	11%	9%	9%	7%	8%

**Agonist/Antagonist: Dorsiflexion/Plantarflexion**

Gender : Age	<30	30-39	40-49	50-59	60-69	>70
male high	27%	27%	33%	30%	28%	29%
male low	28%	29%	29%	28%	31%	39%
average	27%	28%	31%	29%	30%	34%
female high	15%	20%	20%	21%	20%	24%
female low	24%	23%	21%	22%	19%	25%
average	20%	21%	20%	21%	20%	24%

### SHOULDER ABDUCTION/ADDUCTION – 60 degrees per second

BIODEX

\*Reference document: Harbo T, Brincks J, Andersen H (2012) Maximal isokinetic and isometric muscle strength of major muscle groups related to age, body mass, height and sex in 178 healthy subjects. Eur J Appl Physiol 112:267-275

For a complete listing of new normative data.  
Go to [www.biodex.com/clinicalresourcemanual](http://www.biodex.com/clinicalresourcemanual)

# BIODEX

Biodex Medical Systems, Inc.

20 Ramsey Road, Shirley, New York, 11967-4704, Tel: 800-224-6339 (Int'l 631-924-9000), Fax: 631-924-9241, Email: [info@biodex.com](mailto:info@biodex.com), [www.biodex.com](http://www.biodex.com)