

*"The Clinical Advantage"*TM



Biodex Balance System SD

- **Meniscal Repairs**



**Cincinnati SportsMedicine
& Orthopaedic Center**

A Nationally Recognized Center of Excellence

A cooperative effort by Biodex Medical Systems, Inc. and Cincinnati SportsMedicine & Orthopaedic Center

BIODEX

Biodex Medical Systems, Inc.

20 Ramsay Road, Shirley, New York, 11967-4704, Tel: 800-224-6339 (Int'l 631-924-9000), Fax: 631-924-9338, Email: info@biodex.com, www.biodex.com

BALANCE APPLICATION PROTOCOLS

Meniscal Repair

PHASE II – three-four weeks

ACTIVITIES

Weight Shifting.....	4-1
----------------------	-----

PHASE III – five-seven weeks

ACTIVITIES

Weight Shifting.....	4-2
Bilateral Standing / Dynamic.....	4-3

PHASE IV – seven-eight weeks

ACTIVITIES

Bilateral Standing / Dynamic.....	4-4
Single Leg Standing / Static	4-5

PHASE V – nine-twelve weeks

ACTIVITIES

Bilateral Standing / Dynamic.....	4-6
Single Leg Standing / Static	4-7
Postural Stability Test.....	4-7

PHASE VI – thirteen-twenty six weeks

ACTIVITIES

Bilateral Standing / Dynamic.....	4-8
Single Leg Standing / Static	4-9
Postural Stability Test.....	4-9

PHASE VII – twenty seven – fifty two weeks

ACTIVITIES

Bilateral Standing / Dynamic.....	4-10
Single Leg Standing / Static	4-11
Postural Stability Test.....	4-11

Note: The balance progressions that follow are based on the protocols developed and provided by the Cincinnati Sports Medicine Center. The original protocols can be viewed at this link:

<http://www.cincinnati-sportsmed.com/csm/>

All phases are broken down into training and testing possibilities related to these protocols and potential stances utilized.

BIODEX

Biodex Medical Systems, Inc.

20 Ramsay Road, Shirley, New York, 11967-4704

Tel: 800-224-6339 (Int'l 631-924-9000),

Fax: 631-924-9241 Email: info@biodex.com, www.biodex.com

BALANCE APPLICATION PROTOCOLS

Meniscal Repair

PHASE II – THREE-FOUR WEEKS

Balance begins in Phase II: peripheral repairs only

Activities

- Weight Shifting

Weight shifting side to side and forward to back

Use the **Percent Weight Bearing** training screens to perform static medial / lateral weight shifting (*fig.1*), anterior / posterior weight shifting and to re-establish center of balance (*fig.2*).



(fig.1)



(fig.2)

Use **Weight Shift** training screens to emphasize lateral shifting over the affected leg to prepare for full weight bearing ambulation (*fig.3*).



(fig.3)

Positions and Conditions

Bilateral standing / holding on / static mode (*fig.4*).



(fig.4)

BALANCE APPLICATION PROTOCOLS

Meniscal Repair

PHASE III – FIVE-SEVEN WEEKS

Activities

- **Weight Shifting**
- **Bilateral Standing / Dynamic**

Weight shifting side to side and forward to back

Use the **Weight Shifting** training screens in static mode to influence weight bearing to the affected side (*fig.1*) anterior / posterior (*fig.2*) and diagonally to prepare for ambulation (*fig.3*).



(fig.1)



(fig.2)



(fig.3)

Use the **Limits of Stability** training screen in static mode to explore the patient's sway envelope. Train their ability to move away from center, hit a target on the fringe of their allowable sway envelope and return to center (*fig.4*).



(fig.4)

Positions and Conditions

Use bilateral standing / no holding / static mode for **Weight Shifting** training (*fig.5*).

Use bilateral staggered stance / holding on / static mode for diagonal **Weight Shifting** training (*fig.6*).



(fig.5)



(fig.6)

BALANCE APPLICATION PROTOCOLS

Meniscal Repair

Bilateral Standing / Dynamic

Use the **Postural Stability** training screens in dynamic mode to establish postural stability on a moveable surface (*fig.1*) and to have the patient control dynamic movement away from their center of balance (*fig.2*).



(fig.1)



(fig.2)

Positions and Conditions

Use bilateral standing / no holding / dynamic mode for **Postural Stability** training (*fig.3*).



(fig.3)

BALANCE APPLICATION PROTOCOLS

Meniscal Repair

PHASE IV - SEVEN-EIGHT WEEKS

Activities

- Bilateral Standing / Dynamic
- Single Leg Standing / Static

Bilateral Standing / Dynamic

Use the **Postural Stability** training screens in dynamic mode to establish postural stability on a moveable surface (*fig.1*) and to have the patient control dynamic movement away from their center of balance (*fig.2*).



(fig.1)



(fig.2)

Use the **Maze Control** training screen to challenge the patient to control dynamic movement away from their center of balance (*fig.3*).



(fig.3)

Use the **Random Control** training screen to facilitate control of movement around the patient's center of balance which is dictated by the machine (*fig.4*).



(fig.4)

Positions and Conditions

Use bilateral standing / holding on / dynamic mode for Postural Stability, Maze Control and Random Control (*fig.5*).

Progress to bilateral standing / no holding / dynamic mode for Postural Stability, Maze Control and Random Control (*fig.6*).



(fig.5)



(fig.6)

BALANCE APPLICATION PROTOCOLS

Meniscal Repair

Single Leg Standing / Static

Use **Percent Weight Bearing** training screens for single leg activities in static mode to facilitate center of balance on the affected leg medially / laterally (*fig.1*), anteriorly / posteriorly and in combined planes (*fig.2*).



(fig.1)



(fig.2)

Use **Postural Stability** training screens for single leg activity in static mode to facilitate center of balance (*fig.3*).



(fig.3)

Use **Limits of Stability** training screen to challenge the sway envelope of a single leg stance in static mode (*fig.4*).



(fig.4)

Positions and Conditions

Single leg standing / holding (*fig.5*).



(fig.5)

BALANCE APPLICATION PROTOCOLS

Meniscal Repair

PHASE V - NINE-TWELVE WEEKS

Activities

- Bilateral Standing / Dynamic
- Single Leg Standing / Static
- Postural Stability Test (12 weeks)

Bilateral Standing / Dynamic

Use the **Postural Stability** training screens in dynamic mode to establish postural stability on a moveable surface (*fig.1*) and to have the patient control dynamic movement away from their center of balance (*fig.2*).



(fig.1)



(fig.2)

Use the **Maze Control** training screen to challenge the patient to control dynamic movement away from their center of balance (*fig.3*).



(fig.3)

Use the **Random Control** training screen to facilitate control of movement around the patient's center of balance which is dictated by the machine (*fig.4*).



(fig.4)

Positions and Conditions

Use bilateral standing / holding on / dynamic mode for Postural Stability, Maze Control and Random Control (*fig.5*).

Progress to bilateral standing / no holding / dynamic mode for Postural Stability, Maze Control and Random Control (*fig.6*).



(fig.5)



(fig.6)

BALANCE APPLICATION PROTOCOLS

Meniscal Repair

Single Leg Standing / Static

Use **Percent Weight Bearing** training screens for single leg activities in static mode to facilitate center of balance on the affected leg medially / laterally (fig.1), anteriorly / posteriorly and in combined planes (fig.2).



(fig.1)



(fig.2)

Use **Postural Stability** training screens for single leg activity in static mode to facilitate center of balance (fig.3).



(fig.3)

Use **Limits of Stability** training screen to challenge the sway envelope of a single leg stance in static mode (fig.4).



(fig.4)

Positions and Conditions

Single leg standing / holding (fig.5).



(fig.5)

Testing:

Postural Stability / Dynamic Test vs. Normative Data using Fall Risk protocol at 12 weeks (fig.6).



(fig.6)

BALANCE APPLICATION PROTOCOLS

Meniscal Repair

PHASE VI - THIRTEEN-TWENTY SIX WEEKS

Activities

- Bilateral Standing / Dynamic
- Single Leg Standing / Static
- Postural Stability Test (16, 20 and 24 weeks)

Bilateral Standing / Dynamic

Use the **Postural Stability** training screens in dynamic mode to establish postural stability on a moveable surface (*fig.1*) and to have the patient control dynamic movement away from their center of balance (*fig.2*).



(fig.1)



(fig.2)

Use the **Maze Control** training screen to challenge the patient to control dynamic movement away from their center of balance (*fig.3*).



(fig.3)

Use the **Random Control** training screen to facilitate control of movement around the patient's center of balance which is dictated by the machine (*fig.4*).



(fig.4)

Positions and Conditions

Use bilateral standing / holding on / dynamic mode for Postural Stability, Maze Control and Random Control (*fig.5*).

Progress to bilateral standing / no holding / dynamic mode for Postural Stability, Maze Control and Random Control (*fig.6*).



(fig.5)



(fig.6)

BALANCE APPLICATION PROTOCOLS

Meniscal Repair

Single Leg Standing / Static

Use **Percent Weight Bearing** training screens for single leg activities in static mode to facilitate center of balance on the affected leg medially / laterally (fig.1), anteriorly / posteriorly and in combined planes (fig.2).



(fig.1)




(fig.2)

Use **Postural Stability** training screens for single leg activity in static mode to facilitate center of balance (fig.3).



(fig.3)

Use **Limit of Stability** training screen to challenge the sway envelope of a single leg stance in static mode (fig.4).



(fig.4)

Positions and Conditions

Single leg standing / holding (fig.5).



(fig.5)

Testing:

Postural Stability / Dynamic Test vs. Normative Data using Fall Risk protocol at 16, 20 and 24 weeks (fig.6).



(fig.6)

BALANCE APPLICATION PROTOCOLS

Meniscal Repair

PHASE VII - TWENTY SEVEN-FIFTY TWO WEEKS

Activities

- Bilateral Standing / Dynamic
- Single Leg Standing / Static
- Postural Stability Test (28, 32, 36, 40, 44, 48, 52 weeks)

Bilateral Standing / Dynamic

Use the **Postural Stability** training screens in dynamic mode to establish postural stability on a moveable surface (*fig.1*) and to have the patient control dynamic movement away from their center of balance (*fig.2*).



(fig.1)



(fig.2)

Use the **Maze Control** training screen to challenge the patient to control dynamic movement away from their center of balance (*fig.3*).



(fig.3)

Use the **Random Control** training screen to facilitate control of movement around the patient's center of balance which is dictated by the machine (*fig.4*).



(fig.4)

Positions and Conditions

Use bilateral standing / holding on / dynamic mode for Postural Stability, Maze Control and Random Control (*fig.5*).

Progress to bilateral standing / no holding / dynamic mode for Postural Stability, Maze Control and Random Control (*fig.6*).



(fig.5)



(fig.6)

BALANCE APPLICATION PROTOCOLS

Meniscal Repair

Single Leg Standing / Static

Use **Percent Weight Bearing** training screens for single leg activities in static mode to facilitate center of balance on the affected leg medial / laterally (*fig.1*), anterior / posteriorly and in combined planes (*fig.2*).



(fig.1)



(fig.2)

Use **Postural Stability** training screens for single leg activity in static mode to facilitate center of balance (*fig.3*).



(fig.3)

Use **Limits of Stability** training screen to challenge the sway envelope of a single leg stance in static mode (*fig.4*).



(fig.4)

Positions and Conditions

Single leg standing / holding (*fig.5*).



(fig.5)

Testing:

Postural Stability / Dynamic Test vs. Normative Data using Fall Risk protocol at 28, 32, 36, 40, 44, 48, 52 weeks (*fig.6*).



(fig.6)

BIODEX

Biodex Medical Systems, Inc.

20 Ramsay Road, Shirley, New York, 11967-4704, Tel: 800-224-6339 (*Int'l 631-924-9000*), Fax: 631-924-9338, Email: info@biodex.com, www.biodex.com