

WHICH TABLE IS RIGHT FOR YOU?

Urology C-Arm Table - 800

Configure to your exact specifications



SPECIFICATIONS:

Dimensions: 84.5" l x 25.5" w (214.6 x 64.8 cm) with OR accessory rails
 Tabletop: 53" l x 21" w (134.6 x 53.3 cm); 84" l (213.4 cm) with extension
 Lightweight Extension (for patient transfer only, not for imaging): 26.5" l x 21" w (67.3 x 53.3 cm)
 Radiolucent Extension: 31" l x 25.5" w (78.7 x 64.8 cm) with OR accessory rails
 Radiolucent Area: 54" l x 17.5" w (137.2 x 44.5 cm)
 Tabletop Material: Three Fibresin™ sections – each with individual 2" mattress
 Table Sections: Four Sections: Adjustable Fowler back 0° to 80°; abdominal section with cutout; radiolucent removable extension and a lightweight non-imaging transfer extension for transferring patients on/off the table
 Mattress: Individual 2" thick (5 cm), one for each section

Motions:

Height Adjustable: 34" to 43.5" (86.4 to 110.5 cm)
 Tabletop X Motion (head-to-toe): 20" (51 cm)
 Tabletop Y Motion (side-to-side): 10" (25.4 cm)
 Trendelenburg: 0° to 20°
 Reverse Trendelenburg: 0° to 20°
 Lateral Roll: 0° to 20°
 Fowler Back: Adjusts from 0° to 80° with touch handle lock/unlock

Controls:

Hand Control: activates all motions
 Foot Control: activates all motions
 Main Control Panel: LCD screen displays actual tabletop position (all axes) in real time, and provides menus and prompts for other functions
 Joystick Controller: for motorized X-Y positioning

X-Ray Film Tray: Removable film cassette tray for 14" x 17" film with cassette mask for 10" x 12" film. Film cassette tray can be removed from either side of the table.

X-Ray Grid: Removable X-Ray grid, 103 line, 8:1 ratio

Wheels:

Foot End – 3" (7.62 cm) locking swivel casters
 Head End – 5" (12.7 cm) swivel casters with central locking

Accessory Rails: Standard OR accessory rails 1.125" x .375" (2.86 x .95 cm) on table (rails also on table imaging extension)

Patient Restraints: Two body straps, two ankle straps

Finish: Powder coat

Patient Capacity: 400 lb (181.4 kg); weight tested to four times the patient load rating.

Power: 115 VAC or 230 VAC and Battery. Battery automatically charges when table is plugged in. An additional battery can be charged with optional wall-mounted charger.

Warranty: One year

Certifications: UL-2601-1®, CAN/CSA® C22.2 No.: 601-1-M90, EN 60601-1, IEC 60601-1-2:2001-09

058-800 Table, Urology C-Arm, 115 VAC
058-805 Table, Urology C-Arm, 230 VAC

Optional:

056-851 Drain Bag Frame
058-860 Drain Pain, Stainless Steel
056-853 Crutch, Knee, Simplicity, (pair)
058-862 Leg Holder, PAL with boot pad (pair)
058-865 Arm Board, Standard (pair)
058-863 Arm Board, Dura Board with pad (pair)
056-850 IV Pole, rail mounted
056-862 Battery and Wall Charging Unit, 115 VAC
056-863 Battery and Wall Charging Unit, 230 VAC
056-463 Grid, X-ray, 103-line 8:1 ratio, 18" x 18"
058-866 Cassette Tray Holder, 11" x 14"
058-867 Extension Arm Controller

Replacement:

056-857 Drain Bags, Sterile, Disposable, 10/pkg
Fits 056-851 Drain Bag Frame
056-460 Hose, Drain, Urology, 40/cs
056-802 Bellows Skirt, disposable plastic, 250/roll
Fits Surgical C-Arm Tables, 840, Urology C-Arm Table and Brachytherapy Table

Brachytherapy C-Arm Table - 810

Fully image-guided ion implantation



SPECIFICATIONS:

Dimensions: 84.5" l x 25.5" w (214.6 x 64.8 cm) with OR accessory rails
 Tabletop: 53" l x 21" w (134.6 x 53.3 cm); 84" l (213.4 cm) with extension
 Lightweight Extension (For patient transfer only, not for imaging.): 26.5" l x 21" w (67.3 x 53.3 cm)
 Radiolucent Extension: 31" l x 25.5" w (78.7 x 64.8 cm) with OR accessory rails
 Radiolucent Area: 54" l x 17.5" w (137.2 x 44.5 cm)
 Tabletop Material: Three Fibresin™ sections – each with individual 2" mattress
 Table Sections: Four Sections: Adjustable Fowler back 0° to 80°; abdominal section with cutout; radiolucent removable extension and a lightweight non-imaging transfer extension for transferring patients on/off the table
 Mattress: Individual 2" thick (5 cm), one for each section

Motions:

Height Adjustable: 34" to 43.5" (86.4 to 110.5 cm)
 Tabletop X Motion (head-to-toe): 20" (51 cm)
 Tabletop Y Motion (side-to-side): 10" (25.4 cm)
 Trendelenburg: 0° to 20°
 Reverse Trendelenburg: 0° to 20°
 Lateral Roll: 0° to 20°
 Fowler Back: Adjusts from 0° to 80° with touch handle lock/unlock

Controls:

Hand Control: activates all motions
 Foot Control: activates all motions
 Main Control Panel: LCD screen displays actual tabletop position (all axes) in real time, and provides menus and prompts for other functions
 Joystick Controller: for motorized X-Y positioning

X-Ray Film Tray: Removable film cassette tray for 14" x 17" film with cassette mask for 10" x 12" film. Film cassette tray can be removed from either side of the table.

X-Ray Grid: Removable X-Ray grid, 103 line, 8:1 ratio

Wheels:

Foot End – 3" (7.62 cm) locking swivel casters
 Head End – 5" (12.7 cm) swivel casters with central locking

Accessory Rails: Standard OR accessory rails 1.125" x .375" (2.86 x .95 cm) on table (rails also on table imaging extension)

Patient Restraints: Two body straps, two ankle straps

Finish: Powder coat

Patient Capacity: 400 lb (181.4 kg); weight tested to four times the patient load rating.

Power: 115 VAC or 230 VAC and Battery. Battery automatically charges when table is plugged in. An additional battery can be charged with optional wall-mounted charger.

Warranty: One year

Certifications: UL-2601-1®, CAN/CSA® C22.2 No.: 601-1-M90, EN 60601-1, IEC 60601-1-2:2001-09

058-810 Table, Brachytherapy, 115 VAC
058-815 Table, Brachytherapy, 230 VAC

Optional:

056-853 Crutch, Knee, Simplicity, (pair)
058-862 Leg Holder, PAL with boot pad (pair)
058-865 Arm Board, Standard (pair)
058-863 Arm Board, Dura Board with pad (pair)
056-850 IV Pole, rail mounted
056-862 Battery and Wall Charging Unit, 115 VAC
056-863 Battery and Wall Charging Unit, 230 VAC
058-867 Extension Arm Controller

Replacement:

056-802 Bellows Skirt, disposable plastic, 250/roll
Fits Surgical C-Arm Tables, 840, Urology C-Arm Table and Brachytherapy Table