


Which C-Arm Table is right for you?

Surgical C-Arm Table - 840

Three-axis motorized positioning plus X-Y free-float



SPECIFICATIONS:

Dimensions: 104" l x 26" w (260 x 66 cm) with OR accessory rails
 Tabletop: 101" l x 24" w (252 x 61 cm)
 Catheter Tray Extension (removable): 18" l x 24" w (45.7 x 61 cm) non-radiolucent
 Radiolucent Area:
 51" l x 24" w (128 x 61 cm) unobstructed head to accessory rail
 20" l x 17.5" w (50.8 x 44 cm) additional, between accessory rails
 Radiolucent Area, Total Length: 71" l (178 cm)
 Tabletop Material: Carbon fiber with integral head section
 Mattress: 2" thick (5 cm)
 Attenuation:
 Tabletop: .7 mm Aluminum equivalence
 Mattress: .5 mm Aluminum equivalence

Motions:

Height Adjustable: 33.5" to 43.5" (84 to 109 cm)
 Tabletop X Motion (head-to-toe float): 20" (51 cm)
 Tabletop Y Motion (side-to-side float): 10" (25.4 cm)
 Trendelenburg: 0° to 20°
 Reverse Trendelenburg: 0° to 20°
 Lateral Roll: 0° to 20°

Controls:

Hand Control: height, lateral roll and Trendelenburg motions
 Foot Control: activates all motions (has release for float)
 Main Control Panel: LCD screen displays actual tabletop position (all axes) in real time, and provides menus and prompts for other functions
 Push Button: releases X-Y Tabletop motions

Wheels:

Head End - 3" (7.62 cm) swivel casters, integral locking system on base
 Foot End - 5" (12.7 cm) swivel casters with central locking
 Accessory Rails: Standard OR accessory rails 49.5" x 1.12" x .375" (124 x 2.86 x .95 cm) mounted near foot end of table
 Patient Restraints: Three body straps
 Finish: Powder coat
 Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.
 Lift Capacity: May vary between 460 lb - 500 lb (209 - 227 kg) based on line voltage
 Power: 115 VAC or 230 VAC and Battery. Battery automatically charges when table is plugged in. An additional battery can be charged with optional wall-mounted charger.
 Standard Accessory:
 - Catheter Tray Extension (removable)
 - Bellows Skirt (disposable, 250/roll) and dispenser
 Warranty: one year
 Certifications: ETL and cETL Listed to UL-2601-1®, CAN/CSA® C22.2 No.: 601-1-M90, EN 60601-1 and IEC 60601-1-2:2001-09

058-840 Table, Surgical C-Arm - 840, 115 VAC\$33,900.00
058-845 Table, Surgical C-Arm - 840, 230 VAC\$33,900.00

Optional:

056-862 Battery and Charging Unit, 115 VAC.....\$550.00
056-863 Battery and Charging Unit, 230 VAC.....\$550.00
056-850 IV Pole, Rail Mounted\$390.00
056-865 Arm Board, Carbon Fiber, one arm\$600.00
056-866 Arm Board, Carbon Fiber, two arms.....\$1,150.00

Replacement:

056-802 Bellows Skirt, disposable plastic, 250/roll\$95.00
 Fits Surgical C-Arm Tables 830, 840,
 Urology C-Arm Table and Brachytherapy Table

Surgical C-Arm Table - 846

Affordably priced float-top table



SPECIFICATIONS:

Dimensions: 104" l x 26" w (260 x 66 cm) with OR accessory rails
 Tabletop: 101" l x 24" w (252 x 61 cm)
 Radiolucent Area:
 51" l x 24" w (128 x 61 cm) unobstructed head to accessory rail
 20" l x 17.5" w (50.8 x 44 cm) additional, between accessory rails
 Radiolucent Area, Total Length: 71" l (178 cm)
 Tabletop Material: Carbon Fiber with integral head section
 Mattress: 2" thick (5 cm)
 Attenuation:
 Tabletop: .7 mm Aluminum equivalence
 Mattress: .5 mm Aluminum equivalence

Motions:

Height Adjustable: 33.5" to 43.5" (84 to 109 cm)
 Tabletop X Motion (head-to-toe float): 20" (51 cm)
 Tabletop Y Motion (side-to-side float): 10" (25.4 cm)

Controls:

Hand Control: height
 Push Button: releases X-Y Tabletop motions

Wheels:

Head End - 3" (7.62 cm) swivel casters, integral locking system on base
 Foot End - 5" (12.7 cm) swivel casters with central locking
 Accessory Rails: Standard OR accessory rails 49.5" x 1.12" x .375" (124 x 2.86 x .95 cm) mounted near foot end of table
 Patient Restraints: Three body straps
 Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.
 Lift Capacity: May vary between 460 lb - 500 lb (209 - 227 kg) based on line voltage
 Power: 115 VAC or 230 VAC
 Warranty: one year
 Certifications: ETL and cETL Listed to UL-2601-1®, CAN/CSA® C22.2 No.: 601-1-M90 and EN 60601-1

058-846 Table, Surgical C-Arm - 846, 115 VAC\$18,900.00
058-847 Table, Surgical C-Arm - 846, 230 VAC\$18,900.00

Optional:

058-848 Catheter Tray Extension (removable).....\$400.00
056-850 IV Pole, Rail Mounted\$390.00
056-865 Arm Board, Carbon Fiber, one arm\$600.00
056-866 Arm Board, Carbon Fiber, two arms.....\$1,150.00

Which C-Arm Table is right for you?

Pain Management C-Arm Table - 870

Ideal for pain care applications



SPECIFICATIONS:

Dimensions: 93" l x 26" w (233 x 66 cm) with OR accessory rails
 Tabletop: 90" l x 24" w (228.6 x 61 cm)
 Radiolucent Area:
 60" l x 24" w (152.4 x 61 cm) unobstructed head to accessory rail
 4" l x 17.5" w (10.1 x 44 cm) additional, between accessory rails
 Radiolucent Area, Total Length: 64" l (160 cm)
 Tabletop Material: Carbon fiber with integral head section
 Mattress: 2" thick (5 cm)
 Attenuation:
 Tabletop: .7 mm Aluminum equivalence
 Mattress: .5 mm Aluminum equivalence

Motions:

Height Adjustable: 30" to 40" (75 to 100 cm)
 Trendelenburg: 0° to 20°
 Reverse Trendelenburg: 0° to 20°
 Lateral Roll: 0° to 20°

Controls:

Hand Control: activates height, lateral roll and Trendelenburg motions

Optional:

Foot Control: activates height, lateral roll and Trendelenburg motions

Wheels:

Head End – 3" (7.62 cm) swivel casters, integral locking system on base
 Foot End – 5" (12.7 cm) swivel casters with central locking
 Accessory Rails: Standard OR accessory rails 30" x 1.12" x .375"
 (75 x 2.86 x .95 cm) mounted near foot end of table
 Patient Restraints: Three body straps
 Finish: Powder coat
 Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.
 Lift Capacity: May vary between 460 lb - 500 lb (209 - 227 kg) based on line voltage
 Power: 115 VAC or 230 VAC
 Warranty: one year
 Certifications: ETL and cETL Listed to UL-2601-1®, CAN/CSA® C22.2
 No.: 601-1-M90, EN 60601-1 and IEC 60601-1-2

058-870 Table, Pain Management C-Arm - 870,
 115 VAC\$16,900.00
058-875 Table, Pain Management C-Arm - 870,
 230 VAC16,900.00

Optional:

058-877 Controller, Foot\$895.00
056-850 IV Pole, Rail Mounted390.00
056-865 Arm Board, Carbon Fiber, one arm600.00
056-866 Arm Board, Carbon Fiber, two arms.....1,150.00

Surgical C-Arm Table - 846

Affordably priced float-top table



SPECIFICATIONS:

Dimensions: 104" l x 26" w (260 x 66 cm) with OR accessory rails
 Tabletop: 101" l x 24" w (252 x 61 cm)
 Radiolucent Area:
 51" l x 24" w (128 x 61 cm) unobstructed head to accessory rail
 20" l x 17.5" w (50.8 x 44 cm) additional, between accessory rails
 Radiolucent Area, Total Length: 71" l (178 cm)
 Tabletop Material: Carbon Fiber with integral head section
 Mattress: 2" thick (5 cm)
 Attenuation:
 Tabletop: .7 mm Aluminum equivalence
 Mattress: .5 mm Aluminum equivalence

Motions:

Height: Adjustable from 33.5" to 43.5" (84 to 109 cm)
 Tabletop X Motion (head-to-toe float): 20" (51 cm)
 Tabletop Y Motion (side-to-side float): 10" (25.4 cm)

Controls:

Hand Control: height
 Push Button: releases X-Y Tabletop motions

Wheels:

Head End – 3" (7.62 cm) swivel casters, integral locking system on base
 Foot End – 5" (12.7 cm) swivel casters with central locking
 Accessory Rails: Standard OR accessory rails 49.5" x 1.12" x .375"
 (124 x 2.86 x .95 cm) mounted near foot end of table
 Patient Restraints: Three body straps
 Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.
 Lift Capacity: May vary between 460 lb - 500 lb (209 - 227 kg) based on line voltage
 Power: 115 VAC or 230 VAC
 Warranty: one year
 Certifications: ETL and cETL Listed to UL-2601-1®, CAN/CSA® C22.2
 No.: 601-1-M90 and EN 60601-1

058-846 Table, Surgical C-Arm - 846, 115 VAC\$18,900.00
058-847 Table, Surgical C-Arm - 846, 230 VAC18,900.00

Optional:

058-848 Catheter Tray Extension (removable).....\$400.00
056-850 IV Pole, Rail Mounted390.00
056-865 Arm Board, Carbon Fiber, one arm600.00
056-866 Arm Board, Carbon Fiber, two arms.....1,150.00