

"The Clinical Advantage"™

Biodex Offset Unweighing System

Single point suspension permits functional pelvic rotation and versatility when walking, side-stepping, retro-walking and turning

Quick release clips to get a patient out fast, if necessary

Universal Air Support Vest fits a variety of patients

Two position handle for easy and efficient patient unloading

Accommodates children and adults

Attachment rings adjustable for desired degree of pelvic stabilization

Also accommodates other rehab strategies

Attach single side to load involved side



Removable and adjustable patient handrails

Removable Therapist seats provide a place to sit during assistive exercise and provides a safer environment for therapist interaction (not shown) or see page 20.

BIODEX

Dynamic suspension allows for natural vertical displacement of center of gravity

Real time digital scale displays off-loaded weight in pounds or kilograms

Patented off-loading mechanism maintains constant force
Simply dial in amount to offload

Adjust to fit 8' ceilings

Offset frame design offers improved access to manually assist patient footfall and weight shift

Large locking casters for effortless mobility over ground

BIODEX

Biodex Medical Systems, Inc.

20 Ramsay Road, Shirley, New York, 11967-4704

Tel: 800-224-6339 (In NY Call 631-924-9000), Fax: 631-924-9338

Email: sales@biodex.com, www.biodex.com