

BIODEX Gait Trainer 2

Starting Speed - always starts at 0 mph with 0.1 mph speed increments



Heart Rate Monitoring

High Resolution Color Touch-Screen Display

Display Adjustment

Sturdy functional Geriatric Handrails allows clear open platform access

Safety Stop Button

Safety Lanyard

Low Profile

Forward and Reverse Belt Direction

Exact-Track

Elevation 0 to 15% grade

Instrumented Deck - monitors and records step length, step speed and right-to-left time distribution

● Note: Color Printer with Stand - included

BIODEX

Biodex Medical Systems, Inc.

20 Ramsay Road, Shirley, New York, 11967-4704, Tel: 800-224-6339 (Int'l 631-924-9000), Fax: 631-924-9338, Email: info@biodex.com, www.biodex.com