

**CERTIFICATE OF COMPLIANCE FOR RADIOACTIVE MATERIALS PACKAGES  
D.O.T. SPECIFICATIONS 7A, TYPE A**

**Model #001-790 Zevalin Shipping System  
Page 1 of 2**

**1. Package Identification:**

Common Name: Zevalin Shipping System Syringe w/Needle  
Ground or Air Shipment  
Catalog Name: **Shipping System, Zevalin**  
Model Number: **#001-790 with Unit Dose Pigs 001-788 and 01-789**

**2. Preamble:**

- a. This certificate is issued to satisfy section 178.350 including sections 171.2, 173.403, 173.410, 173.411, 173.412, 173.415, 173.431, 173.465, 173.466 and 178.608 of the Department of Transportation Hazardous Materials Regulations (49 CFR).
- b. The packaging and contents described in item 4 below were found to meet the safety standards set forth in, and tested, according to the more stringent recommendations for general performance testing from:  
  
IATA Dangerous Goods Regulations, 58<sup>th</sup> Edition Sections 5.0.4.3, 10.5, and 10.6.1 thru 10.6.3.5.  
  
49 CFR as listed in item 2a above.
- c. This certificate does not release the consignor from compliance with any requirement of the regulations of the U.S. Department of Transportation or other applicable regulatory agencies, including the government of any country through or into which the package will be transported.

**3. Descriptions of Packaging and Authorized contents, Other Conditions, and References:**

- a. **Outer Container:** Model #001-783-A005  
Medium Density Polyethylene: Can hold up to four (4) Zevalin Pigs  
**Dimension:** 11.75" x 11.75" x 12.5" with high density polyurethane foam cushioning and insert support weldment
- b. **Zevalin Unit Dose Pig:** Model #001-788 (system accepts two (2) pigs)  
Lead and acrylic enclosed in 0.065" polycarbonate shell with 0.09" and 0.36" acrylic inside shielding with single turn twist lock top
- c. **Zevalin Unit Dose Pig:** Model #001-789  
0.065" polycarbonate with 0.5" lead body with single turn twist lock top, with 0.42" top lead and 0.56" lead in bottom, and an integral plastic liner
- d. **Absorbent Material:** Model #001-771  
One piece of 3" x 3" material in the bottom of the pig
- e. **Syringe:** B-D or Monoject plastic 10 ml Syringes with needle or luer cap

**CERTIFICATE OF COMPLIANCE FOR RADIOACTIVE MATERIALS PACKAGES  
D.O.T. SPECIFICATIONS 7A, TYPE A**

**Model #001-790 Zevalin Shipping System  
Page 2 of 2**

**Authorized Contents:**

Normal Form Radioactive Materials in Type A quantities

**Specifications and Restrictions:**


Ground or Air Transportation

**4. This certificate is issued on the basis of the test results and safety analysis of the package design.**

**5. Certification**

This document is to certify that to the best of our knowledge, information and belief, the construction methods, the packaging design, and the materials of construction for this specification package comply with the standards described in items 2a, 2b, 2c herein. Any deviation from these conditions must be approved prior to use.

**Approved:**

  
Signature

Clyde Schlein  
Vice President, Compliance & Regulatory Affairs  
Biodex Medical Systems, Inc.

Date: Jan 4, 2013  
Rev: Nov 11, 2015  
Rev: June 7, 2017