

Your Hospital Name  
Address  
City, State Zip  
(000) 000-0000

CHI - SQUARE REPORT

DATE: Dec 01, 1993

TIME: 09:56

ISOTOPE/ACTIVITY: CS-137 / \_\_\_\_\_

COUNT TIME PER INTERVAL: 10 seconds

DETECTOR: PROBE

---

---

INTERVAL #	GROSS COUNTS
1	159744
2	160213
3	159509
4	159445
5	159594
6	160021
7	159701
8	159253 L
9	159488
10	160512 H
	SUM 1597480

  

PROBABILITY	CHI-SQUARE	
0.99	2.09	
0.95	3.33	
0.90	4.17	<lower limit
0.50	8.34	
0.10	14.68	<upper limit
0.05	16.92	
0.01	21.67	

RESULTS: PASSED CHI-SQUARE = 8.25 DEGREES OF FREEDOM = 9

Technologist: \_\_\_\_\_

Comments: \_\_\_\_\_

S/N 871031 D4.0/M4.0