

Which Table is right for you?

Sound Pro Combination Table



SPECIFICATIONS:

Dimensions:

Overall: 71" l x 30" w (180.3 x 76.2 cm); 35" w (88.9 cm) with side rails
 Cardiac Scanning Cutout: 15.75" l x 9.5" w (40 x 24.1 cm)
 Drop-down cushion releases from either side of the table.
 Sonographer Cutout: 14.75" l x 8" w (37.4 x 20.3 cm) Adjustable sonographer cushion drops down or folds up for maximum access to the patient.
 Drop-Down Leg Section: 12.6" l (32 cm)

Motions:

Height Adjustable: 23" to 39" (58.4 to 99 cm)
 Trendelenburg: 0° to ±15°
 Fowler Positioning: 0° to 80° infinitely adjustable

Controls:

Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions
 Foot Control: Activates height, Trendelenburg motions and Fowler positioning
 Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system

Tabletop: Three primary sections: torso section features sonographer and patient cutouts with Fowler positioning; center section remains fixed and leg section drops down to 40° and 80° for stirrup access

Drop-Down Cardiac Scanning Cushion: Release from either side of the table
Adjustable Sonographer Cushion: drops down or folds up for maximum access to the patient

Mattress: Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria
Upholstery Color: Graphite
Patient Restraints: One body strap
Finish: Powder coat
Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.
Weight: 350 lb (159 kg)
Power: 115 VAC or 230 VAC
Certifications: ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1.M90 and EN 60601-1 standards.
Warranty: Two years parts and labor



058-710 Table, Combination, Sound Pro, 115 VAC
Includes: hand and foot controllers, folding side rails, vascular scanning arm board, retractable stirrups, IV pole and paper dispenser.

058-715 Table, Combination, Sound Pro, 230 VAC
Includes: hand and foot controllers, folding side rails, vascular scanning arm board, retractable stirrups, IV pole and paper dispenser.

Ultra Pro Ultrasound Table



SPECIFICATIONS:

Dimensions:

Overall: 70" l x 30" w (177.8 x 76.2 cm); 35" w (88.9 cm) with optional side rails
 Drop-Down Leg Section: 12.6" l (32 cm)

Motions:

Height Adjustable: 23" to 39" (58.4 to 99 cm)
 Trendelenburg: 0° to ±15°
 Fowler Positioning: 0° to 80° infinitely adjustable

Controls:

Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions
 Foot Control (optional): Activates height, Trendelenburg motions, and Fowler positioning
 Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system

Tabletop: Three primary sections: torso section with Fowler positioning, center section remains fixed and leg section drops down to 40° and 80° for stirrup access


Mattress: Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria
Upholstery Color: Graphite
Patient Restraints: Two body straps
Finish: Powder coat
Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.
Weight: 335 lb (152 kg)
Power: 115 VAC or 230 VAC
Certifications: ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1.M90 and EN 60601-1 standards.
Warranty: Two years parts and labor



058-720 Table, Ultrasound, Ultra Pro, 115 VAC
Includes Hand Controller.

058-725 Table, Ultrasound, Ultra Pro, 230 VAC
Includes Hand Controller.

Econo Ultrasound Table



SPECIFICATIONS:

Dimensions:

Overall: 70" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm) with optional side rails
 Drop-Down Leg Section: 12.6" l (32 cm)

Motions:

Height Adjustable: 23" to 39" (58.4 to 99 cm)
 Fowler Positioning: 0° to 80° infinitely adjustable

Controls:

Hand Control: Activates height and Fowler back
 Wheels: 5" (12.7 cm) Individual locking swivel casters

Tabletop: Three primary sections: torso section with Fowler positioning, center section remains fixed and leg section drops down to 40° and 80° for stirrup access

Mattress: Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria
Upholstery Color: Graphite
Patient Restraints: Two body straps
Finish: Powder coat
Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.
Weight: 335 lb (152 kg)
Power: 115 VAC or 230 VAC
Certifications: ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1.M90 and EN 60601-1 standards.
Warranty: Two years parts and labor



058-726 Table, Ultrasound, Econo, 115 VAC
Includes hand controller.

058-727 Table, Ultrasound, Econo, 230 VAC
Includes hand controller.

All ultrasound tables are compatible with Hoyer Lifts

Options are installed by the customer.

Echo Pro Echocardiography Table



SPECIFICATIONS:

Dimensions:

Overall: 73" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm) with optional side rails

Cardiac Scanning Cutout: 15.75" l x 9.5" w (40 x 24.1 cm)

Drop-down cushion releases from either side of the table.

Sonographer Cutout: 14.75" l x 8" w (37.4 x 20.3 cm) Adjustable sonographer cushion drops down or folds up for maximum access to the patient.

Motions:

Height Adjustable: 23" to 39" (58.4 to 99 cm)

Trendelenburg: 0° to ±15°

Fowler Positioning: 0° to 80° infinitely adjustable

Controls:

Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions

Foot Control (optional): Activates height, Trendelenburg motions, and Fowler positioning

Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system

Econo Echocardiography Table



SPECIFICATIONS:

Dimensions:

Overall: 73" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm) with optional side rails

Cardiac Scanning Cutout: 15.75" l x 9.5" w (40 x 24.1 cm)

Drop-down cushion releases from either side of the table.

Sonographer Cutout: 14.75" l x 8" w (37.4 x 20.3 cm) Adjustable sonographer cushion drops down or folds up for maximum access to the patient.

Motions:

Height Adjustable: 23" to 39" (58.4 to 99 cm)

Fowler Positioning: 0° to 80° infinitely adjustable

Controls:

Hand Control: Activates height and Fowler back

Wheels: 5" (12.7 cm) Individual locking swivel casters

Tabletop: Two primary sections: torso section features sonographer and patient cutouts with Fowler positioning; lower body section remains fixed

Mattress: Torso section 2" (5 cm) thick, lower body section 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria

Upholstery Color: Graphite

Patient Restraints: Two body straps

Finish: Powder coat

Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.

Weight: 335 lb (152 kg)

Power: 115 VAC or 230 VAC

Certifications: ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1.M90 and EN60601-1 standards.

Warranty: Two years parts and labor



058-701 Table, Echocardiography, Econo, 115 VAC
Includes band controller.

058-706 Table, Echocardiography, Econo, 230 VAC
Includes band controller.

Echo/Vasc Pro Vascular Echocardiography Table



SPECIFICATIONS:

Dimensions:

Overall: 73" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm) with optional side rails

Cardiac Scanning Cutout: 15.75" l x 9.5" w (40 x 24.1 cm)

Drop-down cushion releases from either side of the table.

Sonographer Cutout: 14.75" l x 8" w (37.4 x 20.3 cm) Adjustable sonographer cushion drops down or folds up for maximum access to the patient.

Motions:

Height Adjustable: 23" to 39" (58.4 to 99 cm)

Trendelenburg: 0° to 15°

Reverse Trendelenburg: 0° to 30°

Fowler Positioning: 0° to 80° infinitely adjustable

Controls:

Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions

Foot Control (optional): Activates height, Trendelenburg motions, and Fowler positioning

Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system

Tabletop: Two primary sections: torso section features sonographer and patient cutouts with Fowler positioning; lower body section remains fixed

Foot Support: Adjustable every 2", up to 14". Stores beneath the table when not in use.

Mattress: Torso section 2" (5 cm) thick, lower body section 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria

Upholstery Color: Graphite

Patient Restraints: Two body straps

Finish: Powder coat

Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.

Weight: 335 lb (152 kg)

Power: 115 VAC or 230 VAC

Certifications: ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1.M90 and EN60601-1 standards.

Warranty: Two years parts and labor



058-702 Table, Vascular Echocardiography, Echo/Vasc Pro, 115 VAC
Includes band controller and foot support.

058-707 Table, Vascular Echocardiography, Echo/Vasc Pro, 230 VAC
Includes band controller and foot support.

Vasc Pro Vascular Ultrasound Table



SPECIFICATIONS:

Dimensions:

Overall: 73" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm) with optional side rails

Motions:

Height Adjustable: 23" to 39" (58.4 to 99 cm)

Trendelenburg: 0° to 15°

Reverse Trendelenburg: 0° to 30°

Fowler Positioning: 0° to 80° infinitely adjustable

Controls:

Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions

Foot Control (optional): Activates height, Trendelenburg motions, and Fowler positioning

Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system

Tabletop: Two primary sections: torso section with Fowler positioning; lower body section remains fixed

Foot Support: Adjustable every 2", up to 14". Stores beneath the table when not in use.

Mattress: Torso section 2" (5 cm) thick, lower body section 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria

Upholstery Color: Graphite

Patient Restraints: Three body straps

Finish: Powder coat

Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.

Weight: 335 lb (152 kg)

Power: 115 VAC or 230 VAC

Certifications: ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1.M90 and EN60601-1 standards.

Warranty: Two years parts and labor



058-732 Table, Vascular Ultrasound, Vasc Pro™, 115 VAC
Includes band controller and foot support.

058-733 Table, Vascular Ultrasound, Vasc Pro™, 230 VAC
Includes band controller and foot support.



058-700 Echocardiography, Echo Pro™, 115 VAC
Includes band controller.

058-705 Table, Echocardiography, Echo Pro™, 230 VAC
Includes band controller.