

The Clinical Advantage™

Ergonomic Sonography Chair



Modern, Ergonomic Design

Sonography Chair features foot ring for additional support and positionings.

www.biodex.com/sonochair

Visit OSHA's website at www.osha.gov on "Sonography: Positioning Patients and Equipment." From sonography chairs to exam tables, Biodex meets all of the sonographer's needs!


BIODEX

www.biodex.com
1-800-224-6339

Int'l 631-924-9000

Ergonomic Sonography Chair

Modern, **ergonomic design** - Sonography Chair features foot ring for additional support and positioning.



The ergonomically designed sonography chair provides multiple adjustments for seat height, back height, and backward tilt tension, along with an integrated seat slider allowing for a variety of comfortable postures.

The comfortable medical-grade vinyl with Advanced Beauty Guard provides protection against bacteria.

A foot ring is standard for additional support and positioning.

SPECIFICATIONS:

- Dimensions:
 - Overall: 23.5" w x 19.5" depth
 - Height Adjustable: 38-47.9
- Seat:
 - Width: 19.3"
 - Depth: 14.5" -18.1"
 - Height: 22.5" - 31.9"
 - Back: 23.5" w x 17.3" h
- Color: Graphite
- Weight Capacity: 250 lb (113.6 kg)
- Shipping Weight: 45 lb (20.4 kg)

058-704 Chair, Ergonomic, Sonography
Includes Foot Ring

www.biodex.com/sonochair

BIODEX
Biodex Medical Systems, Inc.

20 Ramsey Road, Shirley, New York, 11967-4704, Tel: 800-224-6339 (Int'l 631-924-9000), Fax: 631-924-9241, Email: info@biodex.com, www.biodex.com