


Typical Mobile PET Hot Lab



Qty.	Model #	Description
1	086-332	Atomlab 500Plus Dose Calibrator
1	101-356	Source, Cs-137, 200 µCi
1	063-261	Source, Co-57 Simulated Tc-99m, 5 µCi
1	063-139	Rod Source, Cs-137
1	051-013	Model 14C Survey Meter, with Pancake GM Probe
1	019-018	Dosimeter, Personal Electronic, mR
1	042-433	L-Block Shield, Compact, 1.5" lead
1	001-180	Syringe Carrier, Large, .25" lead
1	007-973	Syringe Shield, Pro-Tec PET with lead glass window, 3 cc
1	007-975	Syringe Shield, Pro-Tec PET with lead glass window, 5 cc
1	214-210	Chair, Injection/Resting
1	039-412	Sharps Container Shield, PET, 1" lead
2	039-413	Sharps Container, 3.2 qt., 30/pkg fits 039-412
1	005-400	Radiacwash, Spray Mist, 1 L Bottle
1	121-190	Emergency Spill Kit, Minor
1	024-999	Sign, Caution, Radioactive Materials
1	024-900	Sign, Caution, Radiation Area
1	024-910	Sign, Caution, If You Are Pregnant
1	033-013	Absorbent Paper, Sheets, 50/pkg
2	066-535	Forceps, Straight
1	042-522	Radiation Shield, Mobile, Compact, 1" lead

Options

Qty.	Model #	Description
1	086-338	Dose Calibrator Shielding Rings, Interlocking, 2.25" lead
1	042-434	Lead Brick Cave, 3-wall, 2" lead
1	042-435	Lead Brick Cave Cover
1	002-248	Lead Brick, Rectangular, 8" Long
1	039-106	Waste Container, .125" lead
1	007-980	Syringe Shield, Pro-Tec PET with lead glass window, 10 cc
1	007-945	Syringe Shield, Z-PET, 5 cc
1	086-509	Lineator

Phantoms

Qty.	Model #	Description
1	043-767	Phantom, PET, NEMA 2007/IEC 2008
1	043-768	Phantom, PET Scatter, NEMA
1	043-769	Phantom, PET Sensitivity, NEMA

www.biorex.com/mobilepethotlab

BIODEX

www.biorex.com
1-800-224-6339

Int'l 631-924-9000