

BIODEX

Biodex Medical Systems, Inc.
20 Ramsey Road, Shirley, New York, 11967-4704

FreeStep SAS Supported Ambulation System - Site Survey

Facility and Contact Information

Date:

Facility Name:

Address:

Contact Name:

Phone Number:

Email Address:

Facility Engineer, Maintenance or Operations Person (if any)

Name:

Phone Number:

Email Address:

Site Information

Floor level of installation: Site Access: Elevator: Stairs: Both:

Elevator Type: Maintenance: Passenger:

Elevator Size: L: W: H: Door Opening: W: H:

Doorway or Entrance Way Size: W: H:

Describe receiving capacity. Is there a loading dock? Is there available space to store components until installation? UPS only deliveries?

Is there a requirement for Union labor?

Special Building Codes, regulations, Environmental and/or Hazardous Materials to consider:

Fire Codes: Earthquake: Electrical: Asbestos: Lead Paint:

Description:

Ceiling Information

Drop: Drywall: Open: Other:

Vertical Distance from Floor to Ceiling: Distance from Ceiling to Structural Attachment surface:

Access Description

Describe access to area above the ceiling where the Track will be installed:

Any obstacles?: (Insulation, Oxygen, Water Lines, Ducting, etc...)

Structural Attachment Surface Description

Bar Joist with space: Steel I-Beam: Wood (2" x ____): On Center:

Concrete: Rating:

Services Below Ceiling: *It is best to try to position the **Track** so as not to require moving any of the services below. If they must be moved, it will be the facilities responsibility to do so prior to the Track installation.*

Lights: Fire Suppression: HVAC: Privacy Curtains: Anything Else?:

Description:

Room and Track Layout sketch on the back

Room and Track Layout

Please provide a sketch of the room where the Track is to be installed. Provide relative position to any doors and windows, as well as existing rehabilitation equipment. Include dimensions where possible.

Number of patient support systems:

Each patient support system consists of one (1) trolley, one (1) height adjuster and one (1) spreader bar.
 Note: for 800 Lbs rating - no more than 2 systems per 8 ft of track

Number of harnesses:

At least one per patient support. Additional harnesses are recommended in order to help with patient throughput.

FOR BIODEX USE ONLY					
Support required above ceiling:	Unitrust <input style="width: 50px;" type="text"/>	Threaded Rod: <input style="width: 50px;" type="text"/>	Lateral Bracing: <input style="width: 50px;" type="text"/>		
Total Length: <input style="width: 100px;" type="text"/>	Turntable(s): <input style="width: 50px;" type="text"/>				
8 ft. Straight Lengths: <input style="width: 50px;" type="text"/>	Curves: <input style="width: 50px;" type="text"/>	45°: <input style="width: 50px;" type="text"/>	90°: <input style="width: 50px;" type="text"/>		
10 ft. Straight Lengths: <input style="width: 50px;" type="text"/>	# Brackets Packs: 3 per 8 ft. <input style="width: 50px;" type="text"/>				

Thank you for your interest in the Biodex FreeStep SAS - Supported Ambulation System

Fax completed form to: Biodex Medical Systems, Inc. at (631) 924-9338 or scan and email to sales@biodex.com