



**MIRION**  
TECHNOLOGIES

Accommodates  
**BARIATRIC** Patients

# ECHO PRO™ Echocardiography Table

*Unobstructed patient access ensures sonographer comfort*



## FEATURES

- Ergonomic design helps prevent musculoskeletal injuries – comfortable shoulder, arm and hand positions when scanning
- Extra-wide top with 500 lb patient capacity, suitable for bariatric patients
- Cardiac drop-down cushion, allows release and return from either side of table
- Industry's largest cardiac scanning cutout – 15.75" l x 9.5" w
- Fowler positioning, 0 to 80 degrees, redesigned to accommodate natural body extension and sitting position without slide or shift
- Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties
- Motorized height, Fowler positioning, Trendelenburg/Reverse Trendelenburg and auto-level adjustment
- Naugahyde, antimicrobial mattress cover with Advanced Beauty Gard provides protection against bacteria
- Central floor-locking system maintains table stability
- Durable hand controller activates table positions
- Allows seated or standing work positions

## Optional

- Foot Controller is ideal for when hands are busy with other functions
- Side rails mount flush to table for unencumbered patient access
- Articulating Scanning Arm Board, adjustable from 0 to 130 degrees
- IV Pole easily attaches for quick convenience
- Paper Dispenser accommodates a roll of hygienic table paper. Includes a cutter strap.
- Headrest drops down 15 degrees to allow for carotid and thyroid studies

## DESCRIPTION

The Echo Pro Echocardiography Table is engineered with the sonographer in mind. The ergonomic design provides a barrier-free imaging experience, resulting in less sonographer injury and better quality images.

The table's features make it easier for the sonographer to access the patient from all sides. The cardiac scanning cushion drops down and releases from either side of the table. The sonographer has open access to the left thorax area for an unobstructed apical approach. Side rails mount flush to the table and fold completely beneath the table when not in use.

Patient comfort is assured, even during lengthy procedures with our Naugahyde antimicrobial mattress cover.

Add options such as an Articulating Scanning Arm Board or Headrest for easier thyroid and carotid access to configure the Echo Pro table to fit your department needs.

**BIODEX**  
Part of Mirion Technologies

## Echo Pro Echocardiography Table



The cardiac cushion is designed to release from either side of the table.

### SPECIFICATIONS

#### Dimensions:

- Overall: 73" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm) with optional side rails
- Weight: 335 lb (152 kg)
- Cardiac Scanning Cutout: 15.75" l x 9.5" w (40 x 24.1 cm) drop-down cushion releases from either side of the table
- Sonographer Cutout: 14.75" l x 8" w (37.4 x 20.3 cm) adjustable sonographer cushion drops down or folds up for maximum access to the patient

#### Motions:

- Height Adjustable: 23" to 39" (58.4 to 99 cm)
- Trendelenburg: 0° to +15°, 0° to -15°
- Fowler Positioning: 0° to 80° infinitely adjustable

#### Controls:

- Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions
- Foot Control (optional): Activates height, Trendelenburg motions, and Fowler positioning
- Tabletop: Two primary sections: torso section features sonographer and patient cutouts with Fowler positioning; lower body section remains fixed
- Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system
- Patient Restraints: Two body straps
- Mattress: Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; Naugahyde, antimicrobial mattress cover with Advanced Beauty Gard provides protection against bacteria
- Upholstery Color: Graphite
- Finish: Powder coat
- Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating
- Power: 115 VAC or 230 VAC
- Warranty: Two years parts and labor
- Certifications: ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 6011.M90 and EN 60601-1 standards
- All ultrasound tables are compatible with Hoyer Lifts
- Primary Function: Echocardiography
- CDRH registered and FDA Device Listed



BIO-356 - 11/2021

Copyright © 2021 Mirion Technologies, Inc. or its affiliates. All rights reserved. Mirion, the Mirion logo, and other trade names of Mirion products listed herein are registered trademarks or trademarks of Mirion Technologies, Inc. or its affiliates in the United States and other countries. Third party trademarks mentioned are the property of their respective owners.



Optional flush-mounted side rails provide added safety without obstructing patient access.

### ORDERING INFORMATION

058-700 Table, Echocardiography, Echo Pro, 115 VAC  
Includes hand controller.

058-705 Table, Echocardiography, Echo Pro, 230 VAC  
Includes hand controller.

#### Related:

058-704 Chair, Ergonomic, Sonography

#### Optional:

058-741 Foot Controller

058-633 Side Rails, Folding

058-736 Arm Board, Articulating Scanning

058-737 IV Pole

058-611 Dispenser, Paper

058-738 Headrest

**BIODEX**  
Part of Mirion Technologies

[www.biodex.com](http://www.biodex.com)  
[www.mirion.com](http://www.mirion.com)