

# Biodex Medical Center

55 red st  
Shirley, NY 11967

## Accuracy Test Detail Report

Printed on: 05/08/2017 1:03:56 PM

### Test Information

Test Date	: 03/10/2017 12:20:49 PM	Sealed Source	: Cs-137
Detector S/N	: 96276013	Sealed Source ID	: BM06S-37-6-2
Staff Member	: C S	Calculated Standard	: 174.3 $\mu$ Ci
		Acceptance Variance	: +/- 5 %

### Test Details

Reading	Current Activity	Calculated Activity	Variance
First	176.2 $\mu$ Ci	174.3 $\mu$ Ci	1.1 %
Second	176.6 $\mu$ Ci	174.3 $\mu$ Ci	1.3 %
Third	176.1 $\mu$ Ci	174.3 $\mu$ Ci	1.0 %
Average	176.3 $\mu$ Ci	174.3 $\mu$ Ci	1.1 %

Comments :

Signature : \_\_\_\_\_

Date Signed : \_\_\_\_\_