

# Biodex Medical Center

55 red st  
Shirley, NY 11967

## Geometry Test Detail Report

Printed on: 05/08/2017 1:13:22 PM

### Test Information

Test Date : 03/10/2017 12:23:44 PM      Isotope Name : Tc-99m  
Detector S/N : 96276013      Normalized Volume : 1 ml  
Staff Member : C S

Overall Status: Pass

### Test Details

Container Description:      Type : Vial      Size : 30ml      Material : Plastic

Volume	Measured Activity	Calculated Activity	Correction Factor (Not Applied) 0.95 - 1.05 allowable
1 ml	382. µCi	382. µCi	1.00
2 ml	383. µCi	382. µCi	1.00
5 ml	379. µCi	382. µCi	1.01
10 ml	383. µCi	381. µCi	0.99
15 ml	383. µCi	381. µCi	0.99
20 ml	380. µCi	381. µCi	1.00
25 ml	381. µCi	381. µCi	1.00

Comments :

Signature : \_\_\_\_\_

Date Signed : \_\_\_\_\_