

# BIODEX

Biodex Medical Systems, Inc.  
20 Ramsey Road, Shirley, New York, 11967-4704

## FOR IMMEDIATE RELEASE

### **New Biodex Surgical C-Arm Table 840 Offers Precise Positioning for Image-Guided Fluoroscopic Procedures**

SHIRLEY, N.Y., November 3, 2016 – Biodex Medical Systems, Inc., a provider of innovative medical devices, today announced the availability of the Surgical C-Arm Table 840.

As minimally invasive surgical options have increased over the past decade – with approximately *9 million such procedures* performed in the United States in 2013 – surgeons have come to rely on C-arm fluoroscopy devices for a variety of cardiovascular, orthopedic and urology procedures. The new Surgical C-Arm Table 840 from Biodex is a stable, accessible and vibration-free fluoroscopic table that keeps a patient at exactly the right position and angle for his cardiovascular procedure.

The low-profile base of this radiolucent imaging table from Biodex guarantees C-arm clearance, whether using portable or ceiling-suspended C-arms. It features an extra-large radiolucent area of 71 inches and accommodates patients weighing up to 500 pounds. Proper positioning along three axes is possible with a free-float SmoothGlide™ tabletop movement of 35 inches head to toe and 10 inches side to side. High-speed actuators assure quick tabletop positioning, including height adjustability from 33.5 to 43.5 inches, an isocentric lateral roll of 0° to ±15°, and Trendelenburg motions (and reverse motions) of 0° to 20°. Multiple positioning control access is made possible with mushroom-top, hand and foot controls.

The Biodex Surgical C-Arm Table 840 is available with either a contoured or rectangular carbon-fiber tabletop, both with a 2-inch mattress. The contoured tabletop features a face cutout for comfortable prone positioning. A removable, non-radiolucent catheter tray extension provides additional table length to hold catheters within the patient's sterile field. In order to accommodate a variety of settings, the fluoroscopic table comes standard with both AC and battery power options.

The Biodex Surgical C-Arm Table 840 is part of the Biodex family of C-arm tables. To inquire about this or any model of C-arm table or accessories, please contact Biodex at 800-224-6339 or by email at [info@biodex.com](mailto:info@biodex.com).

#### About Biodex Medical Systems, Inc.

Biodex Medical Systems, Inc. manufactures and distributes medical devices for physical medicine, nuclear medicine, and molecular and medical imaging. It also provides balance assessment for concussion management as well as fall risk screening and conditioning programs. The Biodex commitment to innovative excellence spans more than 60 years. Our customer-driven support is why leading medical facilities around the globe call Biodex first.

###

Contact:

Lila Corwin

VP, Marketing Communications

Biodex Medical Systems, Inc.

20 Ramsey Road

Shirley, NY 11967

631-924-9000, ext. 2281

[lcorwin@biodex.com](mailto:lcorwin@biodex.com)

[www.biodex.com](http://www.biodex.com)