


Biodex Medical Center

55 red st
Shirley, NY 11967

Background Detail Report

Printed on: 05/08/2017 12:55:40 PM



Test Information

Calibration Date/Time : 05/08/2017 10:36:55 AM

Detector S/N : 2007002

Count Time (sec) : 60

Staff : C S

Full Spectrum Information

Total Counts : 470

Counts > 800 keV : 470.0

Total CPM : 470.0

Isotope	ROI (keV)	Counts	CPM	DPM	Net Activity
Wide Window	25 - 800	335	335.0	951	0.00043 μ Ci
Ba-133 wipe	305 - 515	103	103.0	289	0.00013 μ Ci
Ba-wipe-2	256 - 515	123	123.0	282	0.00013 μ Ci
Co-57	50 - 158	79	79.0	90	0.00004 μ Ci
Cs-137	561 - 761	69	69.0	549	0.00025 μ Ci
F-18	420 - 600	74	74.0	271	0.00012 μ Ci
I-131	320 - 410	53	53.0	191	0.00009 μ Ci
Tc-99m	110 - 163	40	40.0	51	0.00002 μ Ci

Comments :

Signature : _____

Date Signed : _____