

# Biodex Medical Center

55 red st  
Shirley, NY 11967

## Chi-Square Detail Report

Printed on: 05/08/2017 12:09:27 PM

### Test Information

Calibration Date/Time : 05/08/2017 10:51:21 AM  
Sample Count Time (sec) : 10

Detector S/N : 2007002  
Staff : C S

### Chi-Square Test Information

Chi-Square Value : 5.43

Status : Pass

Probability	Value
0.99	2.099
0.95	3.33
<b>0.90</b>	<b>4.17</b>
0.50	8.34
<b>0.10</b>	<b>14.68</b>
0.05	16.92
0.01	21.67

(Lower Limit)

(Upper Limit)

Sample	Counts
1	11,715
2	11,760
3	11,696
4	11,812
5	11,916
6	11,786
7	11,695
8	11,821
9	11,836
10	11,934

Comments :

Signature : \_\_\_\_\_

Date Signed : \_\_\_\_\_